

NEW !


iXion® Aligner Pliers

iXion Aligner Pliers are specially designed for effective clear aligner therapy without heating or any aligner distortion. These easy-to-use pliers allow greater torque, rotational control and increased retention.

Use for retention in retainer and aligner systems.


The IX887 Torque Aligner Plier

- Creates a horizontal indentation for individual root
- Indent measures 0.75mm in length and 1.5mm wide




The IX888 Rotation Aligner Plier

- Creates a vertical indentation for individual tooth rotations
- Indent measures 1.5mm in length and 0.75mm wide



The IX889 Retention Aligner Plier

- Creates a circular indentation
- Indent diameter of 1mm



Dystrybucja w Polsce :

Orto Trading, ul. Husarii 11 B, 02-951 Warszawa
tel. 22 651 84 09, 22 497 83 77, 604 234 784